


Product Datasheet

SIVA1 antibody (orb351861)

Description	SIVA1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Apoptosis regulatory protein Siva protein (1-110AA)
Target	SIVA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15304
Dilution Range	IHC-P:1:500-1000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...