


Product Datasheet

TEKT5 antibody (orb351957)

Description	TEKT5 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Tektin-5 protein (1-485AA)
Target	TEKT5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96M29
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH/3T3(Lane 1)...