

Product Datasheet

Omega-gliadin antibody (FITC) (orb352039)

Description	Omega-gliadin antibody (FITC)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	FITC
Immunogen	Recombinant Triticum monococcum Omega-gliadin protein (1-28AA)
Target	Omega-gliadin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02865
Expiration Date	12 months from date of receipt.