



Product Datasheet

ENPP7 antibody (orb35216)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to ENPP7. The protein is a part of the nucleotide pyrophosphatase/phospho... |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | FC, WB |
| Immunogen | Synthetic Peptide |
| Target | ENPP7 |
| Form/Appearance | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Isotype | Rabbit IgG |
| Clonality | Polyclonal |
| MW | 51478 |
| Uniprot ID | Q6UWV6 |
| NCBI | NP_848638.3 |
| Dilution Range | WB: 1:1000, FACS: 1:10-50 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of in HL-60 cell l...