

Product Datasheet

ORF1 antibody (Biotin) (orb352172)

Description	ORF1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Hepatitis E virus genotype 1 Non-structural polyprotein pORF1 protein (60-240AA)
Target	ORF1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P33424
Expiration Date	12 months from date of receipt.