


Product Datasheet

GSTK1 antibody (orb352213)

| | |
|----------------------------|---|
| Description | GSTK1 antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant Human Glutathione S-transferase kappa 1 protein (2-226AA) |
| Target | GSTK1 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q9Y2Q3 |
| Dilution Range | WB:1:1000-5000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of 293T(lane 1), H...