

Product Datasheet

Hemagglutinin antibody (Biotin) (orb352416)

Description	Hemagglutinin antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Rotavirus A Outer capsid protein VP4 protein (247-479AA)
Target	Outer capsid protein VP4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11200
Expiration Date	12 months from date of receipt.