

Product Datasheet

Outer capsid glycoprotein VP7 antibody (orb352453)

Description	Outer capsid glycoprotein VP7 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Rotavirus A Outer capsid glycoprotein VP7 protein (51-326AA)
Target	Outer capsid glycoprotein VP7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10501
Expiration Date	12 months from date of receipt.