


Product Datasheet

IL1R2 antibody (orb352493)

Description	IL1R2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Interleukin-1 receptor type 2 protein (20-320AA)
Target	IL1R2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P27930
Dilution Range	WB:1:200-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse kidney ti...