



Product Datasheet

FKBP4 antibody (orb352513)


Description	FKBP4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase FKBP4 protein (220-459AA)
Target	FKBP4
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q02790
Dilution Range	WB:1:500-2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of using FKBP4 ant...



Immunohistochemical staining of human ep...



Immunoprecipitation analysis of MCF-7 wh...