



Product Datasheet

SIRT3 antibody (orb352737)


Description	SIRT3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human NAD-dependent protein deacetylase sirtuin-3, mitochondrial protein (20-130AA)
Target	SIRT3
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NTG7
Dilution Range	WB:1:200-1000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 whole cel...



Western blot analysis of hepG2 cell usin...



Immunohistochemical staining of human ki...