

Product Datasheet

AMY2B antibody (orb352800)

| | |
|----------------------------|---|
| Description | AMY2B antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC, WB |
| Immunogen | Recombinant Human Alpha-amylase 2B protein (16-320AA) |
| Target | AMY2B |
| Preservatives | PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P19961 |
| Dilution Range | WB:1:200-1000, IHC-P:1:20-1:200 |
| Expiration Date | 12 months from date of receipt. |