


Product Datasheet

PRKG1 antibody (orb352976)

| | |
|----------------------------|---|
| Description | PRKG1 antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant Human cGMP-dependent protein kinase 1 protein (1-300AA) |
| Target | PRKG1 |
| Preservatives | PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q13976 |
| Dilution Range | WB:1:500-2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of mouse lung tiss...