



Product Datasheet

ADCK3 antibody (orb353315)

| | |
|----------------------------|---|
| Description | ADCK3 antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC |
| Immunogen | Recombinant Human Atypical kinase ADCK3, mitochondrial protein (1-220AA) |
| Target | COQ8A |
| Preservatives | PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q8N160 |
| Dilution Range | IHC-P:1:20-1:200 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical staining of human re...



Immunohistochemical staining of human pa...