


Product Datasheet

GLUD2 antibody (orb353452)

Description	GLUD2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Glutamate dehydrogenase 2, mitochondrial protein (319-558AA)
Target	GLUD2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49448
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain tis...