


Product Datasheet

BIRC8 antibody (orb353517)

Description	BIRC8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Baculoviral IAP repeat-containing protein 8 protein (1-236AA)
Target	BIRC8
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96P09
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain tis...