

Product Datasheet

ASCC1 antibody (Biotin) (orb354000)

Description	ASCC1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Activating signal cointegrator 1 complex subunit 1 protein (1-315AA)
Target	Pan myristoylation
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N9N2
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.