

Product Datasheet

Ccr5 antibody (Biotin) (orb354324)

Description	Ccr5 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse C-C chemokine receptor type 5 protein (263-354AA)
Target	Ccr5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51682
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.