


Product Datasheet

CNOT2 antibody (orb354423)

Description	CNOT2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human CCR4-NOT transcription complex subunit 2 protein (1-335AA)
Target	CNOT2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NZN8
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T(Lane 1), J...