

Product Datasheet

FYCO1 antibody (FITC) (orb354985)

Description	FYCO1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human FYVE and coiled-coil domain-containing protein 1 protein (2-255AA)
Target	FYCO1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BQS8
Expiration Date	12 months from date of receipt.