

Product Datasheet

GK5 antibody (HRP) (orb355054)

| | |
|----------------------------|-----------------------------------------------------------------|
| Description | GK5 antibody (HRP) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | HRP |
| Tested Applications | ELISA |
| Immunogen | Recombinant Human Putative glycerol kinase 5 protein (96-301AA) |
| Target | GK5 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q6ZS86 |
| Expiration Date | 12 months from date of receipt. |