



Product Datasheet

HIF1AN antibody (orb355183)


Description	HIF1AN antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Hypoxia-inducible factor 1-alpha inhibitor protein (2-254AA)
Target	HIF1AN
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NWT6
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse heart(Lan...



Immunofluorescence analysis of HeLa cell...



Western blot analysis of MCF-7 whole cel...