

Product Datasheet

HMGCLL1 antibody (orb355199)

Description	HMGCLL1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human 3-hydroxymethyl-3-methylglutaryl-CoA lyase, cytoplasmic protein (2-340AA)
Target	HMGCLL1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TB92
Expiration Date	12 months from date of receipt.