

Product Datasheet

MS4A7 antibody (Biotin) (orb355736)

Description	MS4A7 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Membrane-spanning 4-domains subfamily A member 7 protein (1-47AA)
Target	MS4A7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9GZW8
Expiration Date	12 months from date of receipt.