

Product Datasheet

NKIRAS2 antibody (Biotin) (orb355840)

Description	NKIRAS2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human NF-kappa-B inhibitor-interacting Ras-like protein 2 protein (1-191AA)
Target	NKIRAS2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NYR9
Expiration Date	12 months from date of receipt.