

Product Datasheet

PRPF3 antibody (orb356227)

| | |
|----------------------------|--|
| Description | PRPF3 antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA |
| Immunogen | Recombinant Human U4/U6 small nuclear ribonucleoprotein Prp3 protein (1-200AA) |
| Target | PRPF3 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | O43395 |
| Dilution Range | ELISA |
| Expiration Date | 12 months from date of receipt. |