



Product Datasheet

SLC25A4 antibody (orb356763)


Description	SLC25A4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human ADP/ATP translocase 1 protein (40-74AA)
Target	SLC25A4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P12235
Dilution Range	IHC-P:1:50-200, IF:1:50-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human br...



Immunofluorescence analysis of PC-3 cell...



Immunohistochemical staining of human sm...