


Product Datasheet

SWAP70 antibody (orb357015)

Description	SWAP70 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Switch-associated protein 70 protein (1-300AA)
Target	SWAP70
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UH65
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH/3T3(lane 1)...