

Product Datasheet

UHRF1BP1L antibody (orb357475)

| | |
|----------------------------|--|
| Description | UHRF1BP1L antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA |
| Immunogen | Recombinant Human UHRF1-binding protein 1-like protein (251-522AA) |
| Target | UHRF1BP1L |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | A0JNW5 |
| Dilution Range | ELISA |
| Expiration Date | 12 months from date of receipt. |