

Product Datasheet

ZMYND19 antibody (FITC) (orb357725)

Description	ZMYND19 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Zinc finger MYND domain-containing protein 19 protein (1-227AA)
Target	ZMYND19
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96E35
Expiration Date	12 months from date of receipt.