


Product Datasheet

NGAL antibody (orb357946)


| | |
|----------------------------|---|
| Description | NGAL antibody |
| Species/Host | Mouse |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC, WB |
| Immunogen | Recombinant Human Neutrophil gelatinase-associated lipocalin protein |
| Target | LCN2 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG2b |
| Clonality | Monoclonal |
| Uniprot ID | P80188 |
| Dilution Range | IHC-P: 1:200 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical staining of human sp...



WB analysis of NGAL transfected 293 cell...



Immunohistochemical staining of human pr...