

Product Datasheet

Rabbit Horse IgG antibody (HRP) (orb357972)

Description	Rabbit Horse IgG antibody (HRP)
Species/Host	Rabbit
Reactivity	Equine
Conjugation	HRP
Tested Applications	ELISA
Immunogen	IgG purified from horse serum
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.