


Product Datasheet

SUMO tag antibody (orb357991)

Description	SUMO tag antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	sumo tag protein
Target	Sumo tag
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB:1:500-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of using SUMO tag ...