


Product Datasheet

c Myc antibody (orb357992)

| | |
|----------------------------|---|
| Description | c Myc antibody |
| Species/Host | Mouse |
| Reactivity | All |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, IP, WB |
| Immunogen | EQKLISEEDL (Myc) synthetic peptide conjugate to KLH |
| Target | Myc tag |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Dilution Range | WB:1:500-1:5000, IF:1:100-1:300, IP:1:1000-1:1500 |
| Expiration Date | 12 months from date of receipt. |



Immunoprecipitation analysis of HEK-293

...