


Product Datasheet

Human MAOM protein (Active) (orb359231)

Description	Recombinant human MAOM active protein
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered PBS, pH 7.4.
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His protein
Protein Sequence	MLHIKEKGKPLMLNPRTNKGMAFTLQERQMLGLQGLLPPKIETQDIQALRFHRNL KKMTSPLEKYIYIMGIQERNEKLFYRILQDDIESLMPIVYTPTVGLACSQYGHIFRRP KGLFISISDRGHVRSIVDNWPENHVKAVVVTGGERILGLGDLGVYGMGIPVGKLC LYTACAGIRPDRCLPVCIDVGTDNIALLKDPFYMGLYQKRDRRTQQYDDLIDEFMK AITDRYGRNTLIQFEDFGNHNAFRFLRKYREKYCTFNDDIQGTAVALAGLLAAQ KVISKPISEHKILFLGAGEAALGIANLIVMSMVENGLSEQAQKKIWMFDKYGLLV KGRKAKIDSYPFTHSAPESIPDTFEDAVNILKPSTIIGVAGAGRLFTPDVIRAMA SINERPVI FALS NPTAQAECTAEEAYLTEGRCLFASGSPFGPVKLT DGRVFTPGQ GNNVYIFPGVALAVILCNTRHISDSVFLEAAKALTSQLTDEELAQGRLYPPLANIQE VSINIAIKVTEYLYANKMAFRYPEPEDKAKYVKERTWRSEYDSSL PDVYEWPEAS SPPVITE+HHHHHH
Purity	> 95% as determined by SDS-PAGE and HPLC.
Source	E.Coli
MW	64.4 kDa
Uniprot ID	P23368
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human MAOM protein ...