


Product Datasheet

Human SPRC protein (Active) (orb359252)

Description	Recombinant human SPRC active protein
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered PBS, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His protein
Protein Sequence	MSYYHHHHHHHDYDIPTTENLYFQGAMGSA+PQQEALPDETEVVEETVAEVT EVSVGANPVQVEVGEFDDGAEETEEVVVAENPCQNHHCCKHGKVCELDENN TPMCVCQDPTSCPAPIGEFEKVCSDNKTDFDSSCHFFATKCTLEGTKKGGH KLHLDYIGPCKYIPPCLDSELTEFPLMRDWLKNVLTLYERDEDNLLT EKQKLRVKKIHENEKRLEAGDHPVELLARDFEKNYNMYIFPVHWQFGQLD QHPIDGYLSHTELAPLRAPLIPMEHCTTRFFETCDLDNDKYIALDEWAGC FGIKQKDIDKDLVI
Purity	> 95% as determined by SDS-PAGE and HPLC.
Source	E.Coli
MW	36.1 kDa
Uniprot ID	P09486
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human SPRC protein ...