


Product Datasheet

Human LEG3 protein (Active) (orb359268)

Description	Recombinant human LEG3 active protein
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered 1 × PBS, pH 7.4, 3mM DTT
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Protein Sequence	ADNFSLHDALSGSGNPNPQGWPGAWGNQPAGAGGYPGASYPGAYPGQAPPG AYPGQAPPGAYPGAPGAYPGAPAGVYVPGPPSGPGAYPSSGQPSATGAYPATGP YGAPAGPLIVPYNLPLPGGVVPRMLITILGTVKPNANRIALDFQRGNDVAFHFNPR FNENRRRIVCNTKLDNNWGREERQSVFPFESGKPFKIQVLVEPDHFKVAVNDA HLLQYNHRVKKLNEISKLGISGDIDLTSASYTMI
Purity	> 98% as determined by SDS-PAGE and HPLC.
Source	E.Coli
MW	26.0 kDa
Uniprot ID	P17931
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human LEG3 protein ...