

Product Datasheet

CD177 antibody (FITC) (orb35998)

Description	CD177 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human CD177 antigen protein (22-321AA)
Target	CD177
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N6Q3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.