

---

## Product Datasheet

### EGFR isoform a variant antibody (FITC) (orb360862)

**Description**

EGFR isoform a variant antibody (FITC)

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

FITC

**Tested**

ICC, IF

**Applications**
**Immunogen**

Synthetic peptide (please contact us for details)

**Target**

EGFR isoform a variant (FITC)

**Form/Appearance**

Liquid: PBS, 2 mM Sodium Azide, 30% glycerol

**Concentration**

- 100 µg (in 200 µl): 0.5 mg/ml

**Storage**

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**Purity**

Polyclonal antibodies are purified by peptide affinity chromatography

**Uniprot ID**
[Q9EQ09](#), [P00533](#), [O70156](#)
**NCBI**
[NP\\_005219.2](#), [NP\\_005219.21](#), [NM\\_005228.31](#), [NM\\_005228.3](#)
**Entrez**
[1956](#)
**Dilution Range**

IF/ICC: 1:30-300

**Expiration Date**

12 months from date of receipt.