



Product Datasheet

MRAP2 antibody (orb36099)


| | |
|----------------------------|---|
| Description | MRAP2 antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, IHC, WB |
| Immunogen | Recombinant Human Melanocortin-2 receptor accessory protein 2 protein (66-205AA) |
| Target | MRAP2 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q96G30 |
| Dilution Range | WB:1:1000-5000, IHC-P:1:20-200, IF:1:50-200 |
| Expiration Date | 12 months from date of receipt. |



Immunofluorescence analysis of U251 cell...



Immunohistochemical staining of human br...



Immunohistochemical staining of human lu...