

## Product Datasheet

### CP-544326 (orb363695)

<b>Description</b>	CP-544326 is a selective prostaglandin EP(2) agonist on intraocular pressure in preclinical models o...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{24}H_{22}N_4O_5S$
<b>MW</b>	478.52028
<b>CAS Number</b>	[752187-80-7]
<b>SMILES</b>	<chem>O=S(C1=CN=CC=C1)(N(CC2=CC(OCC(O)=O)=CC=C2)CC3=CC=C(N4C=CC=N4)C=C3)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.