




Product Datasheet

NOTCH3 antibody (orb36386)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	NOTCH3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	243631
Uniprot ID	Q9UM47
NCBI	NP_000426.2
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of H...



Immunohistochemical analysis of formali...