

## Product Datasheet

### Tarafenacin (D-tartrate) (orb363890)

<b>Description</b>	Tarafenacin(SVT-40776) is a highly selective M3 muscarinic receptor antagonist (Ki= 0.19 nM), ~200 f...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{25}H_{26}F_4N_2O_8$
<b>MW</b>	558.4761
<b>CAS Number</b>	[1159101-48-0]
<b>SMILES</b>	<chem>O=C(O[C@H]1CN2CCC1CC2)N(C3=CC=CC(F)=C3)CC4=CC(F)=C(F)C(F)=C4.O[C@@H]([C@H](C(O)=O)O)C(O)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.