

## Product Datasheet

### GSK1059615 (orb363906)

|                        |   |
|------------------------|---|
| <b>Description</b>     | GSK1059615 is a dual inhibitor of PI3K $\alpha$ / $\beta$ / $\delta$ / $\gamma$ (reversible) and mTOR with IC50 of 0.4 nM/0.6 nM/2... |
| <b>Note</b>            | For research use only   |
| <b>Formula</b>         | C <sub>18</sub> H <sub>11</sub> N <sub>3</sub> O <sub>2</sub> S   |
| <b>MW</b>              | 333.3638  |
| <b>CAS Number</b>      | [958852-01-2]   |
| <b>SMILES</b>          | <chem>O=C(S/C1=CC2=CC=C3N=CC=C(C3=C2)C4=CC=NC=C4)NC1=O</chem>   |
| <b>Expiration Date</b> | 12 months from date of receipt.   |