

## Product Datasheet

# Angiotensin I-Converting Enzyme Substrate (orb364247)

<b>Description</b>	Peptide to Angiotensin I-Converting Enzyme Substrate
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Bz-Phe-Ala-Pro
<b>Formula</b>	$C_{24}H_{27}N_3O_5$
<b>Purity</b>	> 95%
<b>Source</b>	Synthetic
<b>MW</b>	437.49
<b>CAS Number</b>	69677-91-4
<b>Expiration Date</b>	12 months from date of receipt.

