

Product Datasheet

Biotin-Exendin 4 (orb364399)

Description	Peptide to Biotin-Exendin 4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Biotin-His-Gly-Glu-Gly-Thr-Phe-Thr-Ser-Asp-Leu-Ser-Lys-Gln-Met-Glu-Glu-Glu-Ala-Val-Arg-Leu-Phe-Ile-Glu-Trp-Leu-Lys-Asn-Gly-Gly-Pro-Ser-Ser-Gly-Ala-Pro-Pro-Pro-Ser-NH ₂
Formula	C ₁₉₄ H ₂₉₆ N ₅₂ O ₆₂ S ₂
Purity	> 95%
Source	Synthetic
MW	4412.96
Expiration Date	12 months from date of receipt.