

Product Datasheet

Deslorelin Acetate (orb364555)

Description	Peptide to Deslorelin Acetate
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Pyr-His-Trp-Ser-Tyr-D-Trp-Leu-Arg-Pro-NH ₂
Formula	C ₆₄ H ₈₃ N ₁₇ O ₁₂
Purity	> 95%
Source	Synthetic
MW	1282.47
CAS Number	57773-65-6
Expiration Date	12 months from date of receipt.