

## Product Datasheet

### Exendin 4 (1-8) (orb364617)

<b>Description</b>	Peptide to Exendin 4 (1-8)
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	His-Gly-Glu-Gly-Thr-Phe-Thr-Ser-NH <sub>2</sub>
<b>Formula</b>	C <sub>35</sub> H <sub>51</sub> N <sub>11</sub> O <sub>13</sub>
<b>Purity</b>	> 95%
<b>Source</b>	Synthetic
<b>MW</b>	833.86
<b>Expiration Date</b>	12 months from date of receipt.