

Product Datasheet

Exendin 4 (3-39) (orb364618)

Description	Peptide to Exendin 4 (3-39)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Glu-Gly-Thr-Phe-Thr-Ser-Asp-Leu-Ser-Lys-Gln-Met-Glu-Glu-Glu-Ala-Val-Arg-Leu-Phe-Ile-Glu-Trp-Leu-Lys-Asn-Gly-Gly-Pro-Ser-Ser-Gly-Ala-Pro-Pro-Pro-Ser-NH ₂
Formula	C ₁₇₆ H ₂₇₂ N ₄₆ O ₅₈ S
Purity	> 95%
Source	Synthetic
MW	3992.47
CAS Number	196109-31-6
Expiration Date	12 months from date of receipt.