




## Product Datasheet

### Pael-R antibody (orb36607)

|                            |                                                                                                                                                                              |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Pael-R                                                                                                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                       |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                 |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                                                                                    |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                                                                                            |
| <b>Target</b>              | GPR37                                                                                                                                                                        |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles                                       |
| <b>Note</b>                | For research use only                                                                                                                                                        |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                   |
| <b>MW</b>                  | 67114                                                                                                                                                                        |
| <b>Uniprot ID</b>          | <a href="#">O15354</a>                                                                                                                                                       |
| <b>NCBI</b>                | <a href="#">NP_005293.1</a>                                                                                                                                                  |
| <b>Dilution Range</b>      | WB: 1:1000, IHC-P: 1:50-100                                                                                                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                              |



Western blot analysis of 293 cell line l...



Western blot analysis of mouse brain tis...